

PROJECT TITLE

SAHAYAK – An AI-Powered External Memory Companion for Alzheimer's & Dementia Patients

An autonomous, wearable AI agent that continuously Observes the environment.

TEAM DETAILS

Team Name: Percevia

Team Members:

- Yash Goel– Team Leader
- Tanish Aggarwal
- Khushi Sharma
- Palak Gupta

College: Vivekananda Institute Of Professional Studies, VIPS-TC

PROBLEM DOMAIN :- Open innovation (Healthcare & IOT)

THE CORE PROBLEM: THE "CONTEXT GAP" & LOSS OF INDEPENDENCE

Alzheimer's and **Dementia** are not merely about forgetting names; they strip away a person's ability to recall the narrative of life, leading to a complete loss of independence.

Patients suffer from **Episodic Memory Loss**, meaning they forget the context of recent events, which manifests in three specific ways:

1. Loss of Object Permanence (Short-Term Lapses):

- Patients constantly struggle with questions like "Where did I put my glasses?" or "Where are my keys?".
- This leads to anxiety and helplessness in daily tasks.

2. Social Disconnection (Prosopagnosia / Face Blindness):




- They may recognize a face but fail to recall the context: "I know him, but is he my son or my neighbor?".
- This causes social withdrawal and isolation.

3. Action & Safety Uncertainty (Visual Agnosia):

- Action Verification: Forgetting critical health tasks like "Did I already take my medicine?".
- Safety Risks: Inability to interpret objects correctly, such as confusing a bottle of cleaning fluid for a drink.

Limitations of Existing Solutions:

- **Smart Speakers (Alexa/Siri):** They are connected to the internet but lack context. They can answer general facts but cannot answer personal questions like "Where is *my* wallet?".
- **GPS Trackers:** They only indicate **location** (where the user is), failing to explain **activity** (what the user is doing).
- **Reminder Apps:** These are **passive** and fundamentally flawed because they require the patient to remember to input the data first.
- **CCTV Cameras:** They provide **passive recording** without intelligence, meaning they cannot process video in real-time or answer user queries.
- **Complex Interfaces:** Many existing apps are too difficult for elderly users to navigate. They need active, voice-guided assistance.
- **Symptom Focus:** Current tools only provide **alerts** rather than actually **recreating the lost memory**.



Proposed Solution: SAHAYAK (THE EXTERNAL HIPPOCAMPUS)



3.1. High-Level Overview

Sahayak is an autonomous, wearable AI agent designed to function as an "Ex Hippocampus" for patients suffering from Alzheimer's and Dementia. Unlike passive GPS trackers or generic smart speakers, Sahayak actively observes the user's environment to create an **Artificial Episodic Memory**. It continuously logs "events" (who, what, when, where) and retrieves this context via voice commands to restore the user's independence.

3.2. Core Philosophy: The Split-Module Architecture

To balance high-performance AI processing with wearer comfort, the solution utilizes a unique "Split-Design" architecture:

- **Module A: The "Eye" (Nural Clip):** A lightweight (<40g), fabric-covered camera and microphone unit worn on the collar. It is designed to look like an accessory, not a medical device.
- **Module B: The "Brain" (Core Unit):** A processing unit (Raspberry Pi 4/5 + Battery) kept in the pocket or a belt pouch. This ensures heat and weight remain away from the user's body.
- **The Link:** A concealed ribbon cable connects the sensory "Eye" to the processing "Brain".

3.3. Technical Implementation Phases

The development is structured in two distinct phases to ensure stability and scalability:

- **Phase 1: Prototype (Laptop-Based):** utilizes a laptop's CPU/GPU for heavy inference (YOLOv8 + CLIP) with a webcam and microphone to validate the software logic.
- **Phase 2: Wearable (Edge-Based):** Migrates the stack to a **Raspberry Pi 4 Model B (4GB)** with a Pi Camera Module 3 and **Bone Conduction Headphones** for private audio output.

SYSTEM ARCHITECTURE (HARDWARE CENTRIC)

4.1. The "Split-Module" Design Philosophy

To address the challenge of running heavy AI models (like YOLO and CLIP) on a wearable device without causing overheating or discomfort, Sahayak employs a **Split-Module Architecture**. This design physically separates the lightweight sensory components (The "Eye") from the heavy processing and power components (The "Brain").

4.2. Module A: The "Neural Clip" (Sensory Unit)

- **Placement:** Worn as a magnetic clip on the shirt collar or chest area.
- **Function:** Acts as the device's "eyes and ears," capturing real-time environmental data.
- **Key Components:**
 - **Visual Sensor:** **Raspberry Pi Camera Module 3** (Wide Angle) is used to capture high-resolution frames for the vision pipeline.
 - **Audio Input:** A **MEMS Microphone** (or USB Mic) captures natural language queries from the user.
 - **Interaction:** Features a tactile **Wake Button** allowing the user to manually trigger the system.
- **Form Factor:** Designed with a "Pebble" aesthetic and fabric cover, weighing less than 40g, to resemble a fashion accessory rather than a medical device.



4.3. Module B: The "Core Unit" (Processing Unit)

- **Placement:** Concealed in the user's pocket or a belt pouch.
- **Function:** Handles all heavy Edge AI processing, power management, and heat dissipation away from the sensitive chest area.
- **Key Components:**
 - **Compute Engine:** **Raspberry Pi 4 Model B (4GB RAM)** (or Pi 5) serves as the central brain handling offline AI inference.
 - **Power Source:** A **10,000mAh Battery** ensures 12+ hours of continuous operation.
 - **Telemetry:** Includes a **GPS Module** for location tracking and a specialized **Heat Dissipation Case** to manage thermal output.

4.4. Interconnectivity & Output

- **The Link:** A thin, concealed ribbon cable runs under the user's shirt, connecting the "Neural Clip" (Eye) to the "Core Unit" (Brain).
- **Audio Output:** The system utilizes **Bone Conduction Headphones** connected via Bluetooth.
 - *Advantage:* These sit on the cheekbones, keeping the ear canal open to ambient sounds, ensuring the user stays aware of their surroundings while receiving private audio cues

FEATURES:-

1. Contextual Social Anchor (Face + History)

- **Function:** The system identifies people and recalls the *context* of the relationship rather than just the name.
- **Output:** Instead of saying "Person Detected," it whispers, *"This is your son, Rahul. He last visited you on Tuesday"*.
- **Benefit:** Addresses **Prosopagnosia** (face blindness) and reduces social anxiety during interactions.

2. Retroactive Memory Search (The Object Finder)

- **Function:** The device passively logs the location of essential items (Keys, Wallet, Glasses) into an offline database using visual object detection.
- **Output:** When the user asks *"Where are my glasses?"*, the system retrieves the last seen location: *"I saw them 15 minutes ago on the kitchen table"*.
- **Benefit:** Solves the issue of **Short-Term Memory Lapses** and object permanence.

3. The "Reality Scanner" (Visual Question Answering)


- **Function:** Uses OCR (Optical Character Recognition) and Visual Question Answering to interpret the safety of objects.
- **Output:** Can read expiry dates or identify dangerous items, warning the user: *"This is a bottle of bleach. Do not drink it"*.
- **Benefit:** Mitigates **Visual Agnosia** and safety risks in daily life.

4. Action Verification & Logging

Function: The system tracks specific actions to verify if a task was completed .

Output: Answers queries like *"Did I take my medicine?"* with specific confirmations: *"Yes, you took the blue pill at 2 PM"* .

Benefit: Prevents accidental overdose or missed medication due to memory gaps.



5. Anxiety Reduction (Voice Cloning)

Function: Utilizes a Text-to-Speech (TTS) engine trained on a loved one's voice (e.g., a family member) instead of a robotic AI voice.

Benefit: Provides emotional reassurance and reduces the patient's fear of using technology.


6. Privacy-First "Shadow" Architecture

Function: All data processing (Vision, Voice, Logic) happens 100% offline on the Raspberry Pi .

Benefit: Ensures **Data Sovereignty**—personal video feeds and memories are never uploaded to the cloud .



Architecture and workflow diagram



1. Input & Perception

Visual: The Camera continuously scans the environment (video feed).

Audio: The Microphone listens for the "Wake Word" or user commands.

2. AI Processing (The "Brain")

Object Detection: YOLOv8 identifies objects (e.g., "Keys", "Medicine").

Identity Matching: FaceNet recognizes people (e.g., "Son", "Doctor").

Contextual Understanding: CLIP verifies if the object is the *same* specific item seen before (e.g., "My specific red cup").

3. Episodic Memory Creation

Stability Check: The system confirms the object is placed (stationary) and not just moving.

Logging: It saves a structured memory log locally: {Object: "Glasses", Location: "Kitchen Table", Time: "10:00 AM"}.

4. Query Processing

User Question: User asks, "Where are my glasses?"

Speech-to-Text: OpenAI Whisper converts the voice into text.

Search: The system scans the local database for the latest entry of "Glasses."

5. Response & Action

Answer Generation: The system formulates a natural sentence: *"You left them on the Kitchen Table."*

Audio Output: Edge-TTS speaks the answer via **Bone Conduction Headphones**.

Advantages of Sahayak

- Fully on-device AI
- Works offline
- Preserves user privacy
- Human-like memory recall
- Wearable & scalable
- Assistive, not intrusive

Applications

- Alzheimer's patients
- Dementia care
- Elderly assistance
- Rehabilitation centers
- Assisted living homes

Limitations

- Limited by Raspberry Pi computer
- Accuracy depends on environment
- Requires calibration per user

Future Scope

- Persistent memory database
- Bone-conduction audio
- Multi-language support
- Caregiver mobile app
- Cloud sync (optional)
- Emotion-aware responses

